

LELAND-JONES PERIPHERAL VASCULAR CLAMP

slightly angled; atraumatic

20.0cm (7.8")



SKU: 142-105-20

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FB468
KLS Martin	24-413-20
AUTHORS	JONES , LELAND
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FB468
KLS Martin	24-413-20
AUTHORS	JONES , LELAND

SPECIFICATIONS

[angled](#)