

LELAND-JONES PERIPHERAL VASCULAR CLAMP

very strong angled, atraumatic

18.0cm (7.0")



SKU: 142-105-18

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FB466
KLS Martin	24-413-18
AUTHORS	JONES , LELAND
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FB466
KLS Martin	24-413-18
AUTHORS	JONES , LELAND

SPECIFICATIONS

[angled](#)