

## LELAND-JONES PERIPHERAL VASCULAR CLAMP

very strong angled, atraumatic

18.0cm (7.0")



**SKU:** 142-105-18

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>Aesculap</b>	FB466
<b>KLS Martin</b>	24-413-18
<b>AUTHORS</b>	<a href="#">JONES</a> , <a href="#">LELAND</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>Aesculap</b>	FB466
<b>KLS Martin</b>	24-413-18
<b>AUTHORS</b>	<a href="#">JONES</a> , <a href="#">LELAND</a>

## SPECIFICATIONS

[angled](#)