

LERICHE HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

15.0cm (6")



SKU: 040-382-15
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH170
KLS Martin	13-328-15
AUTHORS	LERICHE
Medicon	15.58.15
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH170
KLS Martin	13-328-15

AUTHORS

[LERICHE](#)

Medicon

15.58.15

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)