

## LICHTENBERG NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

**SKU:** 050-366-20

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM146
<b>KLS Martin</b>	20-676-20
<b>AUTHORS</b>	<a href="#">LICHTENBERG</a>
<b>Medicon</b>	11.26.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM146
<b>KLS Martin</b>	20-676-20

**AUTHORS**

[LICHTENBERG](#)

**Medicon**

11.26.20

**SPECIFICATIONS**

[serrated](#), [straight](#)