

LIGHT CABLE CONNECTOR

ACMI/Stryker/Codman

unit side

Projector-side
(light input)



SKU: 540-890-01
 Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors