

LIGHT GUIDE CABLE CONNECTOR

Olympus/ACMI - Standard

light source side (old model)

Projector-side
(light input)



SKU: 540-890-10

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors