

## LIGHTLINE TROCAR SLEEVE PLASTIC Ø 11 MM

maintenance free stopcock, WL 90mm

threaded tube, distal bevelled



**SKU:** BLUO/04-00412-1-SF

**Categories:** [Endoscopy](#), [Blue endo](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)