

LOCKHART-MUMMERY FISTULA PROBE

figure 3

16.5cm (6.6")



SKU: 080-130-03

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	fistula probe
Aesculap	BN023
KLS Martin	16-189-16
AUTHORS	LOCKHART-MUMMERY
Medicon	50.94.33



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	fistula probe
Aesculap	BN023
KLS Martin	16-189-16
AUTHORS	LOCKHART-MUMMERY

Medicon

50.94.33