

LOCKHART-MUMMERY FISTULA PROBE

figure 4

16.5cm (6.6")



SKU: 080-130-04

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	fistula probe
Aesculap	BN024
KLS Martin	16-191-16
AUTHORS	LOCKHART-MUMMERY
Medicon	50.94.34



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	fistula probe
Aesculap	BN024
KLS Martin	16-191-16
AUTHORS	LOCKHART-MUMMERY

Medicon

50.94.34