

LOOP ELECTRODE

angled; 27 / 28.5Fr.

brown; for 30° telescopes



SKU: 630-240-01

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes