

LOOP ELECTRODE

angled; 24 / 26Fr.

blue; for 12° telescopes; monopolar



SKU: 630-240-05

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|----------------------------------|
| UNIT | pcs. |
| RANGE | Resectoscopy 4.0 |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|----------------------------------|
| UNIT | pcs. |
| RANGE | Resectoscopy 4.0 |
| TYPE | Electrodes |