

LOUTE WIRE TIGHTENER

21.0cm (8.4")

SKU: 190-780-21

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Wire Twisting Instruments |
| Aesculap | KY220 |
| KLS Martin | 22-542-21 |
| AUTHORS | LOUTE |
| Medicon | 76.05.50 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Wire Twisting Instruments |
| Aesculap | KY220 |
| KLS Martin | 22-542-21 |
| AUTHORS | LOUTE |
| Medicon | 76.05.50 |