

## LOUTE WIRE TIGHTENER

21.0cm (8.4")

**SKU:** 190-780-21

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Twisting Instruments</a>
<b>Aesculap</b>	KY220
<b>KLS Martin</b>	22-542-21
<b>AUTHORS</b>	<a href="#">LOUTE</a>
<b>Medicon</b>	76.05.50

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Twisting Instruments</a>
<b>Aesculap</b>	KY220
<b>KLS Martin</b>	22-542-21
<b>AUTHORS</b>	<a href="#">LOUTE</a>
<b>Medicon</b>	76.05.50