

## M.G.HOSPITAL NEEDLE HOLDER

straight; serrated; grooved jaw

17.0cm (6.8")




**SKU:** 050-072-17

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM257
<b>KLS Martin</b>	20-348-17
<b>AUTHORS</b>	<a href="#">M.G.HOSPITAL</a>
<b>Medicon</b>	10.14.18
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM257
<b>KLS Martin</b>	20-348-17

**AUTHORS**

[M.G.HOSPITAL](#)

**Medicon**

10.14.18

**SPECIFICATIONS**

[serrated](#), [straight](#)