

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

20.0cm (8")



SKU: 050-080-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM237
KLS Martin	20-390-20
AUTHORS	HEGAR , MAYO
Medicon	10.18.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM237
KLS Martin	20-390-20

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.20

SPECIFICATIONS

[serrated](#), [straight](#)