

## MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

20.0cm (8")



**SKU:** 050-080-20

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM237
<b>KLS Martin</b>	20-390-20
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">MAYO</a>
<b>Medicon</b>	10.18.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM237
<b>KLS Martin</b>	20-390-20

**AUTHORS**

[HEGAR](#), [MAYO](#)

**Medicon**

10.18.20

**SPECIFICATIONS**

[serrated](#), [straight](#)