

MC GIVNEY HAEMORRHOIDAL LIGATOR

with 100 ligature rubber rings

Lenght of shaft: 20.0 (8")



SKU: 168-410-00

Categories: [Open Surgery](#), [Stomach, Intestines, Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Haemorrhoidal Ligators
AUTHORS	MC GIVNEY
Medicon	50.93.65

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Haemorrhoidal Ligators
AUTHORS	MC GIVNEY
Medicon	50.93.65