

## MC INTOSH LARYNGOSCOPE SET

3 blades in case; fig. 1-3

fiber optic



**SKU:** 103-300-00

**Categories:** [Anaesthesia](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Fibre Optic Laryngoscopes</a>
<b>Aesculap</b>	AN570
<b>KLS Martin</b>	19-510-00
<b>AUTHORS</b>	<a href="#">MC INTOSH</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Fibre Optic Laryngoscopes</a>
<b>Aesculap</b>	AN570
<b>KLS Martin</b>	19-510-00
<b>AUTHORS</b>	<a href="#">MC INTOSH</a>