

## MC KENZIE SILVER CLIP CARRIER

with clip holding rack

**SKU:** 244-170-01

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Post Mortem Saw</a>
<b>Aesculap</b>	FF421
<b>KLS Martin</b>	24-670-01
<b>AUTHORS</b>	<a href="#">MC KENZIE</a>
<b>Medicon</b>	57.80.11

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Post Mortem Saw</a>
<b>Aesculap</b>	FF421
<b>KLS Martin</b>	24-670-01
<b>AUTHORS</b>	<a href="#">MC KENZIE</a>
<b>Medicon</b>	57.80.11