

## MC KENZIE SILVER CLIP CARRIER

with clip holding rack

**SKU:** 244-170-01

**Categories:** [Open Surgery](#), [Neurosurgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                 |
|-------------------|---------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>            |
| <b>TYPE</b>       | <a href="#">Post Mortem Saw</a> |
| <b>Aesculap</b>   | FF421                           |
| <b>KLS Martin</b> | 24-670-01                       |
| <b>AUTHORS</b>    | <a href="#">MC KENZIE</a>       |
| <b>Medicon</b>    | 57.80.11                        |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                                 |
|-------------------|---------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>            |
| <b>TYPE</b>       | <a href="#">Post Mortem Saw</a> |
| <b>Aesculap</b>   | FF421                           |
| <b>KLS Martin</b> | 24-670-01                       |
| <b>AUTHORS</b>    | <a href="#">MC KENZIE</a>       |
| <b>Medicon</b>    | 57.80.11                        |