

## MENGHINI PUNCTURE NEEDLE

Ø 1.8 x 168mm; Luer-Lock connection

Minimum order quantity: 10 pieces

**SKU:** 086-117-18

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>KLS Martin</b>	18-710-18
<b>Medicon</b>	09.75.78

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>KLS Martin</b>	18-710-18
<b>Medicon</b>	09.75.78