

MENGHINI PUNCTURE NEEDLE

Ø 1.2 x 70mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-03

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF908
KLS Martin	18-700-20
Medicon	09.75.20

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF908
KLS Martin	18-700-20
Medicon	09.75.20