

MENGHINI PUNCTURE NEEDLE

Ø 1.4 x 70mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-05

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF911
KLS Martin	18-700-30
Medicon	09.75.30

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF911
KLS Martin	18-700-30
Medicon	09.75.30