

MENGHINI PUNCTURE NEEDLE

Ø 1.4 x 100mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-06

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF913
KLS Martin	18-700-35
Medicon	09.75.35

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF913
KLS Martin	18-700-35
Medicon	09.75.35