

MENGHINI PUNCTURE NEEDLE

Ø 1.6 x 70mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-07

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF920
KLS Martin	18-700-40
Medicon	09.75.40

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF920
KLS Martin	18-700-40
Medicon	09.75.40