

MENGHINI PUNCTURE NEEDLE

Ø 1.8 x 100mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-08

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Biopsy Needles |
| Aesculap | SF920 |
| KLS Martin | 18-700-50 |
| Medicon | 09.75.55 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Biopsy Needles |
| Aesculap | SF920 |
| KLS Martin | 18-700-50 |
| Medicon | 09.75.55 |