

MENGHINI PUNCTURE NEEDLE

Ø 1.0 x 35mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-01

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF904
KLS Martin	18-700-10
Medicon	09.75.05

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF904
KLS Martin	18-700-10
Medicon	09.75.05