

METAL SHAFT FOR BIPOLAR COAGULATION

forcps 628-101-00 and 628-105-00

Ø 5.0mm; 330mm

No CE mark! No labeling!



SKU: 628-115-53  
 Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

|       |                                         |
|-------|-----------------------------------------|
| UNIT  | <a href="#">pcs.</a>                    |
| RANGE | <a href="#">HF-Instruments, bipolar</a> |
| TYPE  | <a href="#">Sheaths</a>                 |

ZUSÄTZLICHE INFORMATIONEN

|       |                                         |
|-------|-----------------------------------------|
| UNIT  | <a href="#">pcs.</a>                    |
| RANGE | <a href="#">HF-Instruments, bipolar</a> |
| TYPE  | <a href="#">Sheaths</a>                 |