

MICRO DRESSING FORCEPS

straight; serrated; tips 1.0mm

15.5cm (6.2"); width of jaw 1.0mm

SKU: 020-096-10

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	12-184-10
AUTHORS	GERALD
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	12-184-10
AUTHORS	GERALD
SPECIFICATIONS	serrated , straight