

MICRO GRASPER ELECTRODE; BIPOLAR; INSERT

fenestrated; for Ø 5.0mm shaft

for article 628-301-53 and 628-305-53

SKU: 628-341-53

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Instruments, bipolar](#)
TYPE [Electrodes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Instruments, bipolar](#)
TYPE [Electrodes](#)