

## MICRO NEEDLEHOLDER; ROUND HANDLE

straight; 0.8mm; w/o catch

21.0cm (8.4"); diamond coated

**SKU:** 096-200-21

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-012-21
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">without catch</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-012-21
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">without catch</a>