

## MILLER LARYNOSCOPE BLADE

fig. 3; 175mm; fiber optic

**SKU:** 103-322-03

**Categories:** [Anaesthesia](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Fibre Optic Laryngoscopes</a>
<b>KLS Martin</b>	19-585-03
<b>AUTHORS</b>	<a href="#">MILLER</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Fibre Optic Laryngoscopes</a>
<b>KLS Martin</b>	19-585-03
<b>AUTHORS</b>	<a href="#">MILLER</a>