


MILLIGAN SELFRETAINING RETRACTOR

sharp; 3X3 teeth



SKU: 070-840-13

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	MILLIGAN
SPECIFICATIONS	sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	MILLIGAN
SPECIFICATIONS	sharp , toothed