

MINI-HOHMANN BONE LEVER

8mm

16.0cm (6.4")


SKU: 190-571-08

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK118
KLS Martin	23-539-08
AUTHORS	MINI-TOWNSEND



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK118
KLS Martin	23-539-08
AUTHORS	MINI-TOWNSEND