

## MINI-HOHMANN BONE LEVER

8mm

16.0cm (6.4")

**SKU:** 190-571-08

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK118
<b>KLS Martin</b>	23-539-08
<b>AUTHORS</b>	<a href="#">MINI-TOWNSEND</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK118
<b>KLS Martin</b>	23-539-08
<b>AUTHORS</b>	<a href="#">MINI-TOWNSEND</a>