

MIXTER HAEMOSTATIC FORCEPS

curved; serrated

15.5cm (6.2")



SKU: 040-461-15
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-473-15
AUTHORS	MIXTER
Medicon	17.07.36
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-473-15
AUTHORS	MIXTER
Medicon	17.07.36

SPECIFICATIONS

[curved](#), [serrated](#)