

MODELL USA DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 3.4 +/-0.2

SKU: 020-006-15
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-15
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-15
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight