

MODELL USA DRESSING FORCEPS

straight; serrated

30.0cm (12"); width of jaw 5.1 +/-0.2

SKU: 020-006-30
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-31
AUTHORS	STILLE
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-31
AUTHORS	STILLE
SPECIFICATIONS	serrated , straight