

## MODELL USA DRESSING FORCEPS

straight; serrated

12.0cm (4.8"); width of jaw 3.2 +/-0.2

**SKU:** 020-006-12

**Categories:** [Open Surgery](#), [Forceps](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-140-12
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-140-12
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>