

MODELL USA DRESSING FORCEPS

straight; serrated

12.0cm (4.8"); width of jaw 3.2 +/-0.2

SKU: 020-006-12

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-12
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-12
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight