

## MODELL WIEN LIGATURE CATCHER

23.5cm (9.4")

**SKU:** 056-532-23

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Knot Tiers, Ligature Catcher, Loop Twisters</a>
<b>Aesculap</b>	BM952
<b>KLS Martin</b>	20-882-23
<b>AUTHORS</b>	<a href="#">MODELL WIEN</a>
<b>Medicon</b>	13.67.24

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Knot Tiers, Ligature Catcher, Loop Twisters</a>
<b>Aesculap</b>	BM952
<b>KLS Martin</b>	20-882-23
<b>AUTHORS</b>	<a href="#">MODELL WIEN</a>
<b>Medicon</b>	13.67.24