


MONOPOLAR CABLE FOR ELECTRODES AND LOOPS

Ø 4mm; 3.0m for equipm. with banana plug



SKU: 582-533-04
Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Monopolar Cables](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Monopolar Cables](#)