


MONOPOLAR CABLE FOR

Martin units; 3m; Umax: 3 KVp
 instrument and unit side Ø 4.0mm



SKU: 287-815-30
Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables