

MOUSE TEETH FORCEPS

opening to both sides; Ø1.2mm; FL 600mm

4 Fr.; rigid; with LuerLock; insert only



SKU: 583-480-12LLA
Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscscopy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscscopy
TYPE	Semi-rigid Instruments