

MOUSE TEETH/IUD FORCEPS

FL 600 mm; semi-rigid; 5 Fr.

double action; with LuerLock;with handle



with Luer Lock

SKU: 583-480-16LL

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments