

MOUSE TEETH/IUD FORCEPS

FL 600 mm; semi-rigid; 5 Fr.

double action; with LuerLock;insert only



SKU: 583-480-16LLA
Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Uretero-Renoscapy](#)
TYPE [Semi-rigid Instruments](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Uretero-Renoscapy](#)
TYPE [Semi-rigid Instruments](#)