

## MOUSE TEETH/IUD FORCEPS

FL 600 mm; semi-rigid; 5 Fr.

double action; with handle



**SKU:** 583-480-16

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscapy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscapy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>