


MÜLLER APPROXIMATOR

Ø 3.0mm; closing force 25-30g



SKU: 026-530-03

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
AUTHORS	MUELLER
Medicon	56.78.03

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
AUTHORS	MUELLER
Medicon	56.78.03