


## MÜLLER APPROXIMATOR

Ø 4.0mm; closing force 25-30g



**SKU:** 026-530-04

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>AUTHORS</b>	<a href="#">MUELLER</a>
<b>Medicon</b>	56.78.04

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>AUTHORS</b>	<a href="#">MUELLER</a>
<b>Medicon</b>	56.78.04